

Symphony



For more convenience and the adjustment to different applications and skin types, the stitching depth is adjustable.

Micropigmentation is a procedure which is expected to be offered by most beauty salons in the future. The device 'Symphony' is a perfect tool for the practitioner at a very reasonable price.

SAFETY is the most important concern. The 'Symphony' device is operated with the patented safety cartridge system. All parts which could possibly be contaminated are integrated into a single use, sterile cartridge. The cartridge is sealed by a membrane. There is no possibility of backflow of color or body fluids into the handpiece.


PRECISION and **POWER** is necessary to offer a wide range of procedures. The possible applications include lip coloring, eyeliner, fine eyebrow hairstrokes as well as medical procedures and little body tattoos.

NOISE and **VIBRATION** is exceptionally low. A completely new concept has been invented. It is patented. Due to this invention it is possible to manufacture a machine with power and precision without increasing the noise and vibration level. For a price that has not been possible before.

SIZE matters. The handpiece including the needle cartridge is only 12cm (4.7 in.) long, 1,5 cm (0.6 in) in diameter and has a weight of 35g (1.2 oz.).

QUALITY of workmanship is guaranteed with this device. The device 'Symphony' is manufactured in Germany. The R&D and the production is certified according to **ISO 9001** and **ISO 13485** (medical standard).

Attachment of Safety Cartridge



ADVANTAGES on a glance:

- Great Price/ Value Ratio
- Patented Safety Cartridge System
- Adjustable Stitching Depth
- Strong and Precise Mechanism
- Reduced Noise and Vibration
- Exceptional small and lightweight aluminum handpiece
 - Made in Germany
 - 2 years warranty

The **Symphony Set** includes:

- Handpiece
- Footpedal
- Power Supply
- 1, 3, 5 point cartridges (1 each)

